

# Motiva Flora<sup>®</sup> SmoothSilk<sup>®</sup> Tissue Expander Order Form

To Order: 800-924-5072 | [customercare@motivausa.com](mailto:customercare@motivausa.com)

Account #:	
Account Name:	

Motiva Sales Rep:	
P.O. # (if applicable):	

## SURGERY INFORMATION:


Surgery Date:		Patient Reference:	
Shipping Speed:			
Surgeon Name:			
Surgeon Location:			

## CONTACT INFORMATION:

First & Last Name:	
Department:	
Phone & Ext.:	
e-mail:	

## Motiva Flora<sup>®</sup> SmoothSilk<sup>®</sup> Tissue Expander Corsé Projection

Low Height  
Motiva Flora<sup>®</sup> SmoothSilk<sup>®</sup> Tissue Expander  
Single Port, MRI Conditional




Catalog Number	Width (cm)	Height (cm)	Projection (cm)	Volume (cc)	Description	Quantity
L-11-260	11	9	6.7	260	Motiva Flora <sup>®</sup> Corsé Projection Low Height SmoothSilk <sup>®</sup> Expander 260 cc	
L-12-345	12	10	7.1	345	Motiva Flora <sup>®</sup> Corsé Projection Low Height SmoothSilk <sup>®</sup> Expander 345 cc	
L-13-440	13	11	7.6	440	Motiva Flora <sup>®</sup> Corsé Projection Low Height SmoothSilk <sup>®</sup> Expander 440 cc	
L-14-570	14	12	8.0	570	Motiva Flora <sup>®</sup> Corsé Projection Low Height SmoothSilk <sup>®</sup> Expander 570 cc	

Medium Height  
Motiva Flora<sup>®</sup> SmoothSilk<sup>®</sup> Tissue Expander  
Single Port, MRI Conditional



Catalog Number	Width (cm)	Height (cm)	Projection (cm)	Volume (cc)	Description	Quantity
M-11-300	11	10	6.7	300	Motiva Flora <sup>®</sup> Corsé Projection Medium Height SmoothSilk <sup>®</sup> Expander 300 cc	
M-12-375	12	11	7.1	375	Motiva Flora <sup>®</sup> Corsé Projection Medium Height SmoothSilk <sup>®</sup> Expander 375 cc	
M-13-490	13	12	7.6	490	Motiva Flora <sup>®</sup> Corsé Projection Medium Height SmoothSilk <sup>®</sup> Expander 490 cc	
M-14-605	14	13	8.0	605	Motiva Flora <sup>®</sup> Corsé Projection Medium Height SmoothSilk <sup>®</sup> Expander 605 cc	
M-15-750	15	14	8.5	750	Motiva Flora <sup>®</sup> Corsé Projection Medium Height SmoothSilk <sup>®</sup> Expander 750 cc	

High Height  
Motiva Flora<sup>®</sup> SmoothSilk<sup>®</sup> Tissue Expander  
Single Port, MRI Conditional



Catalog Number	Width (cm)	Height (cm)	Projection (cm)	Volume (cc)	Description	Quantity
H-11-345	11	11.5	7.1	345	Motiva Flora <sup>®</sup> Corsé Projection High Height SmoothSilk <sup>®</sup> Expander 345 cc	
H-12-440	12	12.5	7.4	440	Motiva Flora <sup>®</sup> Corsé Projection High Height SmoothSilk <sup>®</sup> Expander 440 cc	
H-13-545	13	13.5	7.5	545	Motiva Flora <sup>®</sup> Corsé Projection High Height SmoothSilk <sup>®</sup> Expander 545 cc	
H-14-680	14	14.5	8.3	680	Motiva Flora <sup>®</sup> Corsé Projection High Height SmoothSilk <sup>®</sup> Expander 680 cc	
H-15-825	15	15.5	8.4	825	Motiva Flora <sup>®</sup> Corsé Projection High Height SmoothSilk <sup>®</sup> Expander 825 cc	
H-16-995	16	16.5	8.8	995	Motiva Flora <sup>®</sup> Corsé Projection High Height SmoothSilk <sup>®</sup> Expander 995 cc	

## Accessory Items

Catalog Number	Description	Quantity
MFPL-001	Motiva Flora <sup>®</sup> Port Locator Replacement	

Motiva Flora<sup>®</sup> SmoothSilk<sup>®</sup> Tissue Expanders are intended for temporary (less than six months) subcutaneous or sub muscular implantation to develop surgical flaps and additional tissue coverage required in a wide variety of applications, particularly to aid in reconstruction following mastectomy, the treatment of underdeveloped breasts and treatment of soft tissue deformities.